PL Savaria portable ceiling lift

| - | - | $\overline{}$ | | | - | Ė |
|---|---|---------------|---|---|---|---|
| Р | ĸ | U | U | U | C | ı |

| Product weight | 5 kg (11 lb) |
|---|--|
| Product life expectancy | 10 years with only 1 battery replacement |
| Safe working load (SWL) | 200 kg (440 lb) or 130 kg (286 lb) |
| Lifting speed | 5.5 cm/sec (2.2 in/sec) at 0 kg, 5 cm/sec (2 in/sec) at 130 kg (286 lb), 4.5 cm/sec (1.8 in/sec) at 200 kg (440 lb) |
| Water ingress protection rating of lift | IP21 |
| Noise level | Maximum 54 dBA |
| Medical electrical equipment class | Class 1 |
| Protection class | Type BF |
| Warranty | 2 years |

BATTERY

| Battery type | Lithium ion, 25.2V 2500 mAh |
|------------------|---|
| | Approximately 30 cycles of 61 cm (24 in) at 200 kg (440 lb) |
| Battery capacity | Approximately 70 cycles of 61 cm (24 in) at 130 kg (286 lb) |
| | Approximately 80 cycles of 61 cm (24 in) at 91 kg (200 lb) |
| Battery charging | Full capacity in approximately 2 hours |

CHARGER

| Charger input | 100-240 VAC, 50-60Hz |
|----------------|----------------------|
| Charger output | 29.4 VDC, max 1A |

ENVIRONMENTAL CONDITIONS

| Operating temperature | 5° C to 40° C (41° F to 104° F) |
|-----------------------|------------------------------------|
| Operating humidity | 15 to 93%, non-condensing |
| Storage temperature | -25° C to 70° C (-13° F to 158° F) |
| Storage humidity | 0 to 93%, non-condensing |
| Atmospheric pressure | 70 kPa to 106 kPA |

CONTROLS

| Lift controls | UP and DOWN buttons to raise and lower the lift |
|---|---|
| Hand control | UP and DOWN buttons to raise and lower the lift |
| Hand control operating force | 4-4.5N |
| Water ingress protection rating of hand control | IP67 |
| On/Off power switch | Used to shut off power to the lift and as an emergency stop feature |

CERTIFICATION

| • ISO 10535:2006 (patient lift) | • IEC 60601-1:2005 A1:2012 (medical electrical equipment) |
|--------------------------------------|---|
| • CAN/CSA Z10535.1:15 (patient lift) | • IEC 60601-1-11:2015 (homecare) |

COMPLIANCE

| COMI EIMICE | |
|---|--|
| • CE Marking per 93/42/EEC (medical device) | • 2011/65/EU (ROHS – 100% of components) |
| • 2006/42/EC (machinery directive) | • 2002/96/EC (WEEE) |

WARNING: This equipment is not suitable in the presence of flammable anesthetic mixtures with air or oxygen, or with nitrous oxide.

